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### Standortinformationen

Shenzhen Watmind Medical Co., Ltd.

8th Floor, Building A, No.16-1, Jinhui Road, Jinsha Community, Kengzishubdistrict, Pingshan District, 518118, Shenzhen, China.

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REF LFA0401-25N

### FOR IN VITRO DIAGNOSTIC USE

This instruction for use (IFU) must be read carefully prior to use. Instruction for use must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions for use.

### PACKING SPECIFICATION

25 Tests/ Kit

### INTENDED USE

This kit is used for in vitro qualitative detection of SARS-CoV-2 antigen in human nasopharyngeal swabs, oropharyngeal swabs and saliva samples.

Positive result from the test need further analyze with clinical history of patient and other diagnostic information to determine patient infection status. Positive value is only a reference guide for clinical diagnosis. The test results only reflect the current state of the sample. Negative result cannot exclude SARS-CoV-2 infection and should not be used as the sole basis for patient management decisions. Negative results must be combined with clinical observations, patient history, and epidemiological information.

This test is only for clinical laboratory, medical institutions use or for real-time inspection by professional medical staff, not for home testing. It cannot be used as the basis for the diagnosis and exclusion of pneumonia caused by SARS-CoV-2.

The laboratory testing of SARS-CoV-2 should meet the requirements of the "Laboratory testing for SARS-CoV-2 in suspected human cases" and other requirements, and pay attention to the biosecurity.

### PRINCIPLE OF THE PROCEDURE

The SARS-CoV-2 Ag Diagnostic Test Kit (Colloidal Gold) is a colloidal gold enhanced double antibody sandwich immunoassay for the qualitative determination of Nucleocapsid (N) Protein of SARS-CoV-2 virus in human nasopharyngeal swabs, oropharyngeal swabs and saliva samples.

SARS-CoV-2 antibody is immobilized in the test region on nitrocellulose membrane. If the specimen contains SARS-CoV-2 antigen, during the assay specimen is allowed to react with the colored conjugate (SARS-CoV-2 antibody-colloidal gold conjugate); the mixture then migrates chromatographically on the membrane by the capillary action. An SARS-CoV-2 positive specimen produces a distinct color band in the test region, formed by the specific antibody antigen colored conjugate complex "(Au-SARS-CoV-2-Ab)-(SARS-CoV-2-Ag)-(SARS-CoV-2-Ab)". Absence of this colored band in the test region suggests a negative result. A colored band always appears in the control region serving as procedural control regardless of the specimen contains SARS-CoV-2 or not.

### REAGENTS AND MATERIALS SUPPLIED

#### 1. Main components:

Specification	components
Ingredients	
Test card with desiccant in a sealed foil pouch	25
Extraction tube holder	1
Tube with sample extraction solution	25
Nasopharyngeal swab (with thin head & soft pole)	25
Oropharyngeal swab (with thick head & hard pole)	25
Instruction for use	1

#### 2. Ingredients of the Test device

SARS-CoV-2 antibody	Coated in the Test region on NC membrane
Goat anti Chicken IgY polyclonal antibody	Coated in the control region on NC membrane
SARS-CoV-2 antibody, Chicken IgY, colloidal gold conjugate	Coated in the conjugate pad
Other test device supports	/

#### 3. Ingredients of the sample extraction solution

- Phosphate solution

Note: The components in different batches of the kit cannot be mixed.

### MATERIALS REQUIRED BUT NOT PROVIDED

- Timer
- Pipette

### STORAGE AND STABILITY

- Kits should be stored in 2°C~30°C in a cool, dark, dry place preservation, valid for 18 months (Tentative), forbidden to store under 2°C and avoid using expired products.
- The test card should be in aluminum foil bag after opening, to the specified environment (temperature 2°C~35°C, humidity 40%~60%) used within 15 minutes.
- The buffer should be used immediately after dropping into the dropper.
- MFD date and EXP date: marked on the label.

### SPECIMEN REQUIREMENTS

- Collection of nasopharyngeal secretion: Insert the sterile swab into the place where the nasopharyngeal secretions are the most, and rotate the swab close to the inner wall of the nasal cavity 3 times, remove the swab.
- Collection of oropharyngeal secretion: Insert the sterile swab from mouth completely into the oropharyngeal swelling, centering on the red part of the throat wall, epicondylitis, and tonsils, wipe and rotate 3 times with moderate force to avoid touching the tongue and remove the swab.
- Collection of saliva samples: Press the tip of tongue against the root of jaw to concentrate saliva. Apply the Oropharyngeal swab under the tongue for at least 10 seconds and soak it completely.
- The samples should be used as soon as possible after collected (within half an hour).
1. Nasopharyngeal and oropharyngeal swabs are stable in 30 mins when kept in the sample extraction solution provided with the kit.  
Samples should be tested as soon as possible after collection.
2. If transport of samples with universal transport medium (UTM-RT® System, Copan

Diagnostics, Murrieta, CA, USA) is required, minimal dilution of the sample is recommended, as dilution may result in decreased test sensitivity. Whenever possible, 1 mL or less is best to avoid excessive dilution of the patient sample. Nasopharyngeal and oropharyngeal swabs in UTM are stable for up to 48h at 2°C-8°C.

- Samples should not be inactivated.

### TEST PROCEDURE

Do not open the pouch until you are ready to perform a test, and the single-use test is suggested to be used under low environment humidity (RH≤70%) within 1 hour.

- Allow all kit components and specimens to reach room temperature between 18°C~26°C prior to testing.
- Remove the test card from the foil pouch and place on a clean dry surface.
- Identify the test card for each specimen.


### Sample processing:

- Elute swab with Sample extraction solution

1.1 Place the swab into the sample tube that has been pre-filled with sample extraction solution and then completely immerse the swab head in the sample tube. Vigorously mix the solution by rotating the swab forcefully against the side of the tube at least 10 times (while submerged) and squeeze the tube 5 times by hand to ensure that the sample on the sampling swab is fully eluted into the sample extraction buffer.

1.2 Squeeze the swab head along the inner wall of the extraction tube to keep the liquid in the tube as much as possible. Discard the swab and cover the drip head to mix the liquid thoroughly.

\*Note: Recommend to use a pipette to transfer the samples to reduce deviations.

2. Prepare swabs in UTM (UTM-RT® System, Copan Diagnostics, Murrieta, CA, USA): Remove 2ml or more UTM, 1 mL or less left in the tube is best to avoid excessive dilution of the patient sample. Insert the swab into the tube until the breakpoint is level with the tube opening and rotate the swab at least 10 times. Bend the swab shaft at a 180 degrees angle to break it off at the breaking point. You may need to gently rotate the swab shaft to complete the breakage. Samples are stable for up to 48h at 2°C-8°C.

\*Note: When using UTM, it is important to ensure that the UTM containing the sample is warmed to room temperature (18°C~26°C). Cold samples will not flow correctly and can lead to erroneous or invalid results. Several minutes will be required to bring a cold sample to room temperature.

### Test operation

Before performing the test, you must read the instruction manual of the product completely, and please balance the test cards and sample extraction solution to room temperature (18°C~26°C) before the test. Do not perform the test only when the reagent was equilibrated to room temperature (18°C~26°C) to avoid affecting the accuracy of the experimental results.

3. Remove the test card from the foil pouch and place on a clean dry surface. Dispense 100µL (3drop) of the specimen into the circular sample well on the card.
4. Interpret the test results at 10-15 minutes. Do not interpret the results after 20 minutes.

Discard used test tubes and Test device in suitable biohazards waste container Caution: Use a clean pipette or tip for every sample to avoid cross-contamination.

### INTERPRETATION OF TEST RESULTS

**Positive:** Both purplish test band and purplish control band appear on the membrane.

**Negative:** Only the purplish control band appears on the membrane. The absence of a test band indicates a negative result.

**Invalid:** There should always be a purplish control band in the control region regardless of test result. If control band is not seen, it indicates that the incorrect operation process or the kit has deteriorated or damaged or the antibody content in the specimen is too high. In this case, read the instructions carefully again and dilute the sample to be tested with a new test device. If the problem persists, stop using this lot number immediately and contact your local supplier.

Note: The purplish band in the test area (T) can show the color depth. However, within the specified observation time, regardless of the color of the ribbon, even a very weak ribbon should be judged as a positive result.

### LIMITATIONS

- The result of the product should not be taken as a confirmed diagnosis, for clinical reference only. Judgement should be made along with RT-PCR results, clinical symptoms, epidemic condition and further clinical data.
- Negative results, from patients with symptom onset beyond seven days, should be treated as presumptive and confirmation with a molecular assay, if necessary, for patient management, may be performed.
- A negative test result may occur if the level of antigen in a sample is below the detection limit of the test.

4. Due to the limitation of the detection method, the negative result cannot exclude the possibility of infection. The positive result should not be taken as a confirmed diagnosis. Judgement should be made along with clinical symptoms and further diagnosis methods.
5. This reagent can only qualitatively detect SARS-CoV-2 antigens in human nasopharyngeal swab, oropharyngeal swab and saliva samples. It cannot determine the certain antigen content in the samples.
6. The accuracy of the test depends on the sample collection process. Improper sample collection, improper sample transportation and storage or freezing and thawing of the sample will affect the test results.
7. It is optimum when eluting swabs with the matched samples extraction solution. Using other diluents may result in wrong results.
8. The solution and test card must be equilibrated to room temperature (18°C~26°C) before used, otherwise the results may be incorrect.
9. Sensitivity maybe decrease if the sample did not test directly. Please test the sample as soon as possible.
10. Cross reactions maybe exist due to the N protein in SARS has a high homology with the new coronavirus (SARS-CoV-2). However, the interpretation of the results is not affected during seasons without SARS infection.
11. The following reasons may cause false negative results:
  - 1) Inappropriate sample collection, using other non-matching solution, sample transfer time is too long (more than half an hour), the volume of solution added when eluted the swab are too much, non-standardized elution operation, low virus titer in the sample, these may all lead to false negative results.
  - 2) Mutations in viral genes may lead to changes in antigen epitope, leading to false negative results.
12. Analysis the possibility of false positive results:
  - 1) Inappropriate sample collection, using other non-matching solutions, non-standardized elution operation, these may all lead to false positive results.
  - 2) Cross-contamination of samples may lead to false positive results.
  - 3) False negative result from nucleic acid.
13. Analysis the possibility of invalid result:
  - 1) If the sample volume is not enough, the chromatography cannot be carried out successfully.
  - 2) The test card would invalid if the package was broken. The packaging status must be carefully checked before use.

#### PERFORMANCE CHARACTERISTIC

1. Corporate control
  - 1.1. The coincidence rate of positive controls  
Tested with 5 positive controls, the results were all positive, and the coincidence rate (+ / +) was 5/5.
  - 1.2. The coincidence rate of negative controls  
Tested with 7 negative controls, the results were all negative, and the coincidence rate (+ / +) was 7/7.
- 1.3. Repeatability  
Tested with repeatable control 1 for 10 times, the results were all positive and consistent. Tested with repeatable control 2 for 10 times, the results were all positive and consistent.
- 1.4. Limit of detection  
Use 3 different concentration LoD controls to test, L1 is negative, L2 ~ L3 are positive. Note control samples, L1~L3 are all internal corporate control. This kit was confirmed to detect  $1.5 \times 10^2$  TCID<sub>50</sub>/mL of SARS-CoV-2 which was isolated from USA-WA1/2020, Gamma-Irradiated.
- 1.5. Cross reactivity  
The following viruses and other Microorganism have no effect on the test results.

Potential Cross-Reactant	Test Concentration	Test Result	
Virus	Respiratory Syncytial Virus A	$1.0 \times 10^5$ PFU/mL	No cross reaction
	Respiratory Syncytial Virus B	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No cross reaction
	Measles Virus	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No cross reaction
	Adenovirus Type 3	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No cross reaction
	Adenovirus Type 7	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No cross reaction
	Human cytomegalovirus	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No cross reaction
	Varicella-zoster virus	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No cross reaction
	Human coronavirus OC43	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No cross reaction
	Human coronavirus 229E	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No cross reaction
	Rotavirus	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No cross reaction
	Influenza B	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No cross reaction
	Influenza A	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No cross reaction
	Mycoplasma pneumonia	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No cross reaction
	Epstein Barr Virus	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No cross reaction
	MERS-CoV	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No cross reaction
Coronavirus NL63	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No cross reaction	
Other Microorganism	Staphylococcus aureus	$1.0 \times 10^6$ CFU/mL	No cross reaction
	Streptococcus pneumoniae	$1.0 \times 10^6$ CFU/mL	No cross reaction

- 1.6. Interfering Substances  
The following interfering substances have no effect on the test results.

Substance	Active Ingredient	Concentration	Test Result
Endogenous	Mucin	2.0 % w/v	No interference
Nasal Spray	Oxymetazoline	12 % v/v	No interference
Sore Throat Phenol Spray	Phenol	15 % v/v	No interference
Throat Lozenge	Benzocaine, Menthol	0.15% w/v	No interference
Anti-viral Drug	Ribavirin	12.9 mg/mL	No interference
Antibacterial, Systemic	Tobramycin	4.0 ug/mL	No interference

#### 2. Clinical study

##### 2.1 Nasopharyngeal swabs and oropharyngeal swabs samples

The clinical evaluation was performed to compare the results obtained with the SARS-CoV-2 Ag Diagnostic Test Kit and a comparative reverse transcriptase polymerase chain reaction test (Novel Coronavirus (2019-nCoV) Nucleic Acid Diagnostic Kit (PCR-Fluorescence Probing) manufactured by Sansure Biotech Inc). Patients who were treated within of 14 days showed symptoms were included in the study. All patients come from the same hospital. By RT-PCR 185 positive and 228 negative samples, 224 samples were from nasopharyngeal and 189 samples from oropharyngeal swabs. For the nasopharyngeal swabs samples, 103 were positive and 121 were negative. For the oropharyngeal swabs samples, 82 were positive and 107 were negative. The presentation of the results of the SARS-CoV-2 Ag Diagnostic Test Kit components based on the different parameters is as follows summarized:

CT value	Number of samples	2019 nCoV RT-PCR Results	SARS-CoV-2 antigen test result as compared to RT-PCR
≤30	107	positive	103/107=96.26% (95%CI:90.78%-98.54%)
≤36	185	positive	161/185=87.03% (95%CI:81.42%-91.13%)
>40	228	negative	226/228=99.12% (95%CI:96.86%-99.76%)

days	Number of samples	2019 nCoV RT-PCR Results	SARS-CoV-2 antigen test result as compared to RT-PCR
≤7	103	positive	98/103=95.15% (95%CI:89.14%-97.91%)
≤14	133	positive	122/133=91.73% (95%CI:85.80%-95.32%)
>14	52	positive	39/52=75.00% (95%CI:61.79%-84.77%)

Sensitivity: 87.03% (95%CI:81.42%~91.13%) for CT values ≤36  
Sensitivity: 95.15% (95%CI:89.14%~97.91%) for onset of symptoms within 7 days  
Specificity: 99.12% (95%CI:96.86%~99.76%)

##### 2.2 Saliva samples

The clinical evaluation was performed to compare the results obtained with the SARS-CoV-2 Ag Diagnostic Test Kit and a comparative reverse transcriptase polymerase chain reaction test (Novel Coronavirus (2019-nCoV) Nucleic Acid Diagnostic Kit (PCR-Fluorescence Probing) manufactured by Sansure Biotech Inc). Among patients, there are 157 positive and 182 negative samples by RT-PCR confirmed. The presentation of the results of the SARS-CoV-2 Ag Diagnostic Test Kit is as follows:

CT value	Number of samples	2019 nCoV RT-PCR Results	SARS-CoV-2 antigen test result as compared to RT-PCR
≤30	46	positive	44/46=95.65% (95%CI:85.47%-98.80%)
≤36	157	positive	142/157=90.45% (95%CI:84.84%-94.12%)
>40	182	negative	181/182=99.45% (95%CI:96.95%-99.90%)

days	Number of samples	2019 nCoV RT-PCR Results	SARS-CoV-2 antigen test result as compared to RT-PCR
≤7	89	positive	85/89=95.51% (95%CI:89.01%-98.24%)
≤14	116	positive	104/116=89.66% (95%CI:82.79%-93.98%)
>14	41	positive	35/41=85.37% (95%CI:71.56%-93.12%)

Sensitivity: 90.45% (95%CI:84.84%-94.12%) for CT values ≤36  
Sensitivity: 95.51% (95%CI:89.01%~98.24%) for onset of symptoms within 7 days  
Specificity: 99.45% (95%CI:96.95%~99.90%)

#### PRECAUTIONS

1. The reagent is a disposable diagnostic reagent in vitro, which is only used for the detection of human nasopharyngeal swab, oropharyngeal swab and deep sputum samples. The operation should be carried out strictly according to the instructions. Do not use expired and damaged products.
2. The kit should be sealed and kept away from moisture. Reagents or samples stored at low temperature should be balanced to room temperature before they can be used.
3. Reagents should be used as soon as possible after removal from aluminum foil bags, so as to avoid exposure to air for too long and affecting test results due to dampness.
4. Do not use samples that have been placed for too long or contaminated.
5. Please operate in accordance with the laboratory testing procedures for infectious diseases. Waste after use should be treated in accordance with infectious substances and should not be discarded at will.
6. Incorrect operation may affect the accuracy of the results, such as insufficient sample mixing, insufficient amount, inaccurate detection time, etc.
7. Components in different batch should not be mixed.
8. There should be appropriate biosafety assurance procedures for those substances containing and suspected sources of infection. The following are relevant considerations:
  - 1) Handle samples and reagents with gloves;
  - 2) Do not suck samples with your mouth;
  - 3) Do not smoke, eat, drink, cosmetic or handle contact lenses while handling these items;
  - 4) Disinfect the spilled sample or reagent with disinfectant;
  - 5) Disinfect and treat all samples, reagents and potential pollutants in accordance with relevant local regulations;

6) Each component of the reagent remains stable until the expiry date under proper handling and storage conditions. Do not use the expired reagent kit.

**BIBLIOGRAPHY**

- 1, Peaper DR, Landry ML. Rapid diagnosis of influenza: state of the art. Clin Lab Med. 2014;34(2):365 - 385.
2. Patel J, Sharma P. Design of a novel rapid immunoassay for simultaneous detection of hepatitis C virus core antigen and antibodies. Arch Virol. 2020;165(3):627 - 641.
3. Chafekar A, Fielding BC. MERS-CoV: Understanding the Latest Human Coronavirus Threat. Viruses. 2018;10(2):93. Published 2018 Feb 24.

**[Explanation of Labels]**

 IVD	In Vitro Diagnostic Use	 See Instruction for Use	 REF	Catalog #
 LOT	Batch Number	 Expiry Date		Manufacturing Date
 Do not reuse		 Store between 2~30°C		Keep away from Sunlight
 Keep Dry		 Manufacturer		EU Authorized Representative
 CE Mark				

 Shenzhen Watmind Medical Co., Ltd.

8th Floor, Building A, No.16-1, Jinhui Road, Jinsha Community, Kengzi subdistrict, Pingshan District, 518118, Shenzhen, China.  
 Tel: +86 755-86969964  
 Fax: +86 755-26658059  
 Website: <http://www.watmind.com>

 Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537, Hamburg, Germany  
 Tel: +49-40-2513175  
 Fax: +49-40-255726

**[Date of Approval and Amendment of IFU]: 16 Jan, 2021**